PassportTM/HydroportTM * Single-piece Foldable Lens Placement System

DESCRIPTION
The PassportTM single-piece foldable lens placement system is a device used for folding and delivering the Bausch & Lomb one-piece C11UB lens into the eye. The Hydroport single-piece foldable lens placement system is a device used for folding and delivering the Bausch & Lomb one-piece BL-27[™] or Akreos[®] Adapt[™] lenses into the eye. The devices provide a small tubular pathway in which the C11UB, BL-27 or Akreos Adapt lenses can be placed into the eye with one continuous forward motion. The Passport/Hydroport single-piece foldable lens placement system consists of two parts:

1. Syringe shaped tube (includes lens loading deck and lens cap) and plunger (Fig. 1).
 2. Transition cell (Fig. 2).
-


The Passport/Hydroport single-piece foldable lens placement system is a sterile (EO), disposable plastic device, designed for single use only.

INDICATIONS


The PassportTM single-piece foldable lens placement system should only be used with the Bausch & Lomb C11UB series intraocular lenses. The Hydroport single-piece foldable lens placement system should only be used with the Bausch & Lomb BL-27 and Akreos Adapt series intraocular lenses from 10 to 30 diopters in power.

INSTRUCTIONS FOR USE

1. Open the package using standard sterile procedures. Place contents onto the sterile field.
2. Fill the rectangular end of the transition cell with approximately 0.2 mL of a Bausch & Lomb approved viscoelastic, using the line as a guide (Fig. 3). NOTE: Ocucoat[™] viscoelastic is validated for use with the Hydroport single piece foldable lens placement system.


3. Select the syringe-shaped tube and plunger. Pull the plunger until the first O-ring or gasket is just inside of the syringe-shaped tube. **Verify that the tip of the plunger is even with or behind the loading deck.** Place a drop of viscoelastic directly in the center of the loading deck and in the fork of the plunger. With a blunt, non-toothed forceps, place the lens into the recessed area in the loading deck (Fig. 4a or 4b).

4. Advance the plunger until the fork fully captures the haptic or for the Akreos Adapt model the lens edge. Place a drop of viscoelastic in the center of the lens. Verify that the lens is inside of the recessed area of the loading deck, prior to closing the lens cap.
5. Close the lens cap (Fig. 5). Make sure no part of the lens is pinched between the loading deck and cap. If this is the case, open the lens cap and reposition the lens.
6. Align the key on the lens cap with the slot in the transition cell (Fig. 6). Push together with a firm motion until the tube with the cap and transition cell snap together (Fig. 7).

Viscoelastic material will flow and fill both the lens chamber and the transition cell. **Once assembled, pull on the transition cell to verify that it cannot be separated from the tube.**

7. Fill the distal end of the device with viscoelastic material or balanced salt solution to reduce the possibility of introducing air pockets into the eye during implantation.

8. The lens is now ready to be delivered. Turn the device over, so that the bevel is face down, and insert the beveled portion of the transition cell through the incision in the eye. When the tip is positioned, apply continuous uniform forward motion to the plunger. **A slow smooth forward motion of the plunger is essential for proper transition of the lens into position. Reverse movement of the plunger may result in the loss of the control of the lens or damage to the lens.** Continue forward motion until the lens is fully expressed from the tip.

9. The haptic of the unfolded lens may remain engaged in the plunger fork, or in the case of the Akreos Adapt model the lens edge until it is released into the eye. To disengage the lens, release pressure on the plunger end allowing it to retract into the tip.
10. Discard the device after use.

PRECAUTIONS

Store the Passport/Hydroport single-piece foldable lens placement system's sterile sealed package at or below 30°C/86°F. Contents are sterile unless package is opened or damaged. Do not resterilize. Do not reuse. Single use only. Discontinue lens delivery if excessive force is needed to advance the plunger and lens.

MEDICAL DEVICE RE-USE STATEMENT

If this product is re-processed and/or re-used, Bausch & Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as "single use" which is defined as a device intended to be used once only for a single patient.

SURGICAL PROCEDURE

Proper surgical technique is the responsibility of the individual surgeon. Surgeons must determine the suitability of any particular procedure based upon their medical training and experience.

WARRANTY AND LIMITATIONS OF LIABILITY

Bausch & Lomb warrants that the Passport/Hydroport single-piece foldable lens placement system, when delivered, will conform to the manufacturer's then current version of the published specifications for such Passport/Hydroport single-piece foldable lens placement system in all material respects and shall be free from defects in material or workmanship for a period equal to the Passport/Hydroport single-piece foldable lens placement system's expiration date. **Bausch & Lomb excludes all other warranties, whether expressed, implied or by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for a particular purpose. Bausch & Lomb shall not be liable for any incidental, consequential or exemplary loss, damage or expense, directly or indirectly arising from the use of the Passport/Hydroport single-piece foldable lens placement system even if Bausch & Lomb has been advised of the possibility of such loss, damage or expense.**

RETURN GOODS POLICY

All products returned to Bausch & Lomb Incorporated must be accompanied by a Return Goods Authorization Number. Call USA 1-800-338-2020. Outside the USA call your customer service affiliate.

ADVERSE REACTION REPORTING

Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as Passport/Hydroport related and that were not previously expected in nature, severity or degree of incidence should be reported to Bausch & Lomb at 1-800-356-9382. Outside the USA call your customer service affiliate.

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