BAUSCH + LOMB

enVista[™]TORIC

HYDROPHOBIC ACRYLIC IOL

en


Device Description

The enVista[™] toric hydrophobic acrylic intraocular lens (IOL) (non-preloaded model: ETE / preloaded into shuttle models: ETPE, ETC20E, ETC24E) was developed to replace the natural crystalline lens in adult patients in whom the cataractous lens has been removed. The composition and characteristics of the IOL are specified in the table below.

The enVista toric IOL has an aspheric optic that is designed to be free of spherical aberration through the incorporation of a proprietary spherical aberration-neutral optic design that does not influence the pre-surgical corneal spherical aberration profile. The optic is designed with the SureEdge[™] posterior-squared step edge to provide a 360-degree PCO barrier. The IOL employs an Accuset[™] haptic with a broad, modified C-loop design and optic-haptic offset to facilitate improved contact and stability within the capsular bag. The posterior-located cylinder axis marks denote the meridian with the lowest power. The enVista material that makes up the TruSight[™] optic has been assessed for glistening-free capacity and scratch resistance. Incorporated StableFlex[™] technology allows for enhanced IOL compliance for ease of loading, premium control through IOL delivery, and efficient optical recovery following implantation.

Physical Characteristics

Lens / Haptic Material	Hydrophobic acrylic (hydroxyethyl methacrylate (HEMA)-polyethylene glycol phenyl ether acrylate (poly(EG)PEA)-styrene copolymer, crosslinked with ethylene glycol dimethacrylate)									
Material Characteristics	Index Of Refraction: 1.53 @ 35 °C; Specific Gravity: 1.19 g/ml									
Optic Type / Powers	Aspheric / +6.0 to +30.0 Diopters in 0.5 Diopter increments (SE – Spherical Equivalent)									
Cylinder Powers (CYL) – IOL Plane	0.90 D ¹	1.25 D	1.50 D	2.00 D	2.50 D	3.00 D	3.50 D	4.25 D	5.00 D	5.75 D
Cylinder Powers (CYL) – Corneal Plane	0.64 D	0.90 D	1.06 D	1.40 D	1.76 D	2.11 D	2.45 D	2.98 D	3.50 D	4.03 D
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 12.5 mm; Haptic Angle: 0°									
Spectral Transmittance	Ultraviolet: UV (389) 10 % transmittance for +20.0 Diopter IOL. See Figure 1 with chart's X value = Wavelength (nm) and Y value = % Transmittance; chart compares the transmittance curve of an enVista IOL to a 53-Year-Old Human Lens. A = +20 Diopter enVista IOL and B = 53-Year-Old Human Lens. NOTE: Light transmittance values for an IOL material may vary slightly depending on the method of measurement. Reference: 53-Year-Old Human Lens data from Boettner, E.A. and Wolter, J.R., "Transmission of the Ocular Media," Investigative Ophthalmology, 1:776-783, 1962.									
	¹ Only available in Canada									

Indications

The enVista toric hydrophobic acrylic IOL (non-preloaded model: ETE / preloaded into shuttle models: ETPE, ETC20E, ETC24E) is indicated for primary implantation in the capsular bag of the eye in adult patients for the visual correction of aphakia and corneal astigmatism following the removal of a cataractous lens for improved uncorrected distance vision.

Warnings

As with any surgical procedure, there is risk involved. Physicians considering IOL implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

1. Recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
4. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage to the endothelium during implantation.
6. Suspected microbial infection.
7. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.
8. Rotation of the enVista toric IOL away from the intended axis can reduce its astigmatic correction. Misalignment greater than 30° may increase postoperative refractive cylinder. If necessary, IOL positioning should occur prior to capsule fibrosis and IOL encapsulation.

Precautions For Use And Storage

1. Do not attempt to resterilize the IOL as this can produce undesirable side effects.
2. Prior to opening, inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use. The IOL should be used immediately after opening.
3. Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage (such as the loss of saline storage solution, or the presence of salt crystallization).
4. Store at room temperature. Do not freeze. Avoid high temperatures (>43 °C / >109 °F). Keep dry. Keep away from sunlight.
5. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution or sterile normal saline.
6. Do not place the IOL in contact with surfaces where such contamination can occur.
7. Do not autoclave the IOL.
8. Do not re-use the IOL. It is intended for permanent implantation. If explanted, sterility and proper function cannot be ensured.
9. The safety and effectiveness of the IOL have not been substantiated in patients with pre-existing ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting an IOL in a patient with one or more of these conditions. Physicians considering IOL implantation in such patients should explore the use of alternative methods of aphakic correction and consider IOL implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.

Before Surgery

- Retinal conditions, or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment

or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this IOL

- Amblyopia
- Clinically severe corneal dystrophy (e.g., Fuchs')
- Rubella, congenital, traumatic or complicated cataracts
- Extremely shallow anterior chamber, not due to swollen cataract
- Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye (e.g., iritis or uveitis)
- Aniridia
- Iris neovascularization
- Glaucoma (uncontrolled or controlled with medication)
- Microphthalmos or macrophthalmos
- Optic nerve atrophy
- Previous corneal transplant
- Pre-existing ocular conditions which may negatively impact stability of the implant

During Surgery

- Mechanical or surgical manipulation required to enlarge the pupil
 - Vitreous loss (significant)
 - Anterior chamber bleeding (significant)
 - Uncontrollable positive intraocular pressure
 - Complications in which the IOL stability could be compromised
10. Patients with preoperative problems, such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, and chronic drug miosis may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from IOL implantation when such conditions exist.
 11. A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on IOL implantation before attempting to implant the IOL.
 12. As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block, cyclitic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, acute corneal decompensation, toxic anterior segment syndrome (TASS), and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: IOL repositioning, IOL replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.
 13. Care should be taken to remove all viscoelastic from the anterior and posterior surfaces of the lens to minimize the possibility of lens rotation causing misalignment of the IOL from the intended axis placement.
 14. Effectiveness of implanting a toric IOL in reducing postoperative astigmatism is affected by many factors, including the following:
 - The degree of mismatch between the postoperative magnitude of corneal

astigmatism and effective IOL power in the corneal plane.

- Misalignment between the intended axial position and final IOL axial orientation.
- Error in prediction of the postoperative corneal cylinder axis and power. Error in prediction of cylinder axis is greatest for lower levels of preoperative corneal astigmatism.
- Manufacturing variation in power and axis markings can influence intended correction.

Medical Device Re-Use Statement

If this product is reprocessed and/or re-used, Bausch + Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as 'single-use' which is defined as a device intended to be used once only for a single patient.

Calculation Of IOL Power

Suggested A-Constant: 119.1 (OPTICAL BIOMETRY)

The recommended A-Constant is intended for use with axial length measurements obtained by optical biometry. Use of axial length measurements by other techniques (e.g., Applanation A-scan) will normally require a different IOL constant. This number is a guideline only and is based on an evaluation of clinical data obtained using the IOL Master. The physician should determine preoperatively the power of the IOL to be implanted.

Selection And Placement Of enVista Toric IOL

Keratometry and biometry data should be used in place of refractive data to determine targeted amount of astigmatism correction. Incision size and location influence amount of postoperative corneal astigmatism and its respective axis. It is recommended that surgeons customize their surgically-induced corneal astigmatism values based upon individual surgical technique and past results. To facilitate IOL selection and axis placement, Bausch + Lomb provides a web-based proprietary tool, the enVista toric Calculator (<https://envista.toriccalculator.com>), for the surgeon. Preoperative keratometry and biometry data, incision location, spherical equivalent IOL power, and the surgeon's estimated surgically-induced corneal astigmatism are used as inputs for the enVista toric Calculator. These inputs are used to determine the axis of placement in the eye and the predicted residual refractive astigmatism for up to three different enVista toric IOL models. In eyes with low levels of corneal astigmatism, the predicted residual refractive astigmatism for implantation of the enVista toric IOL will be displayed for evaluation by the surgeon to determine the clinically meaningful benefit of implanting a toric IOL. Prior to surgery, the operative eye should be marked in the following manner: With the patient sitting upright, precisely mark the twelve o'clock and/or the six o'clock position with a T marker, a surgical skin marker, or a marking pencil indicated for ophthalmic use. Using these marks as reference points, an axis marker can be used immediately prior to or during surgery to mark the axis of IOL placement. Input from the enVista toric Calculator can be used to determine optimal axis of placement. Toric axis markings at the haptic-optic junction identify the flat meridian of the enVista toric IOL and represent an imaginary line of the plus cylinder axis. After the IOL is inserted in the capsular bag, precisely align the axis markings on the enVista toric IOL with the marked axis of IOL placement. Be sure to remove all viscoelastic from the capsular bag. Reconfirm proper alignment of the enVista toric IOL following

viscoelastic removal and/or inflation of the capsular bag at the end of the surgical case. Residual viscoelastic and/or over-inflation of the bag may cause lens rotation away from the intended axis of placement. Deviation from the intended axis of placement may compromise effectiveness of astigmatic correction. Inaccurate astigmatism measurements, errors in corneal markings, inaccurate placement of the enVista toric IOL axis during surgery, unanticipated surgically-induced changes in the cornea, or physical rotation of the lens after implantation may also limit the desired effect of the toric IOL on correction of corneal astigmatism.

Directions For Use

1. Inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use.
2. Prior to implanting, examine the IOL package for type, power, and proper configuration.
3. Open the peel pouch and remove the vial in a sterile environment.
4. Remove the lid from the vial.
5. Follow steps below.
 - a) Non-preloaded IOL (model: ETE)
 1. With a pair of smooth forceps, remove the IOL from the vial by gently grasping the IOL haptic.
 2. Rinse the entire IOL with sterile balanced salt solution or sterile normal saline.
 3. Examine the IOL thoroughly to ensure particles have not become attached to it, and examine the lens optical surfaces for other defects.
 4. The IOL may be soaked in sterile balanced salt solution until ready for implantation.
 5. It is recommended to use an approved inserter per the Validated Inserters table below.
 - b) Preloaded IOL (models: ETPE, ETC20E, ETC24E)
 1. Remove the SnapSet™ IOL shuttle from the vial by grasping the tab on the top.
 2. The SnapSet IOL shuttle may be flushed with sterile balanced salt solution during preparation for delivery.
 3. The SnapSet IOL shuttle preloaded with the enVista IOL (models: ETPE, ETC20E, ETC24E) is designed for use with the EyeGility™ inserter for enVista preloaded.
6. It is recommended to use an approved viscoelastic for lubrication of the IOL during implantation. See table below.
7. There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.

Validated Inserters

Model	Inserter	Viscoelastic
ETE	BLIS	Amvisc™ Plus, Amvisc™, OcuCoat™
	INU100	
ETPE (IOL only)	EyeGility (2.0 mm or 2.4 mm)	
ETC20E (IOL + inserter)	EyeGility 2.0 mm	
ETC24E (IOL + inserter)	EyeGility 2.4 mm	

Overview Of Clinical Studies

The enVista model ETE is a result of minor modifications (that did not require clinical study data) from the parent enVista model

MX60 and its predecessor enVista model MX60T.

A clinical study of the MX60 began in the United States on October 19, 2010. This prospective, single arm, open label study included a total of 122 subjects (122 eyes) at 6 clinical sites. Postoperatively, subjects underwent complete ophthalmic evaluations at regularly scheduled intervals through Form 4 (Postoperative Days 120-180). At the Form 4 visit, 118 subjects (100 %) achieved BCVA of 20/40 or better, which exceeds the FDA grid of 96.7 %. The rates of FDA defined potentially sight-threatening adverse events that occurred in the clinical trial at Form 4 were found to be less than the "FDA Grid" of Historical Controls. Two cumulative adverse events (2/122; 1.6 %) of cystoid macular edema were reported through the Form 4 visit. One persistent adverse event (1/121; 0.8 %) of cystoid macular edema was reported at the Form 4 visit. No serious ocular adverse events occurred during this study. All subjects in the safety analysis set were evaluated for IOL glistenings at Form 3 and Form 4 visits. IOL glistenings were evaluated via retroillumination slit lamp examination utilizing a photographic grading scale provided in the protocol. The grading scale consisted of (in order of severity), "none, grade 0 (trace), grade 1, 2, 3, or 4." No glistenings of any grade were reported for any subject at any visit in the clinical study. The results of clinical investigation provided reasonable assurance that the Model MX60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

enVista MX60T Clinical Trial Study Results

The US clinical trial of the enVista toric intraocular lens was conducted in 191 subjects (191 eyes). The dioptric power range was 16.0 to 27.0 D with cylindrical powers at the lens plane of 1.25 D, 2.00 D, and 2.75 D for the MX60T.

Study Description

The study was a prospective, multicenter, parallel-group, partially randomized, partially controlled, double-masked, monocular clinical trial to evaluate the safety and effectiveness of the enVista toric IOL, Model MX60T, in reducing postoperative refractive astigmatism in subjects undergoing cataract extraction. Subjects in the lowest astigmatic IOL power (1.25 D) cohort were randomized to undergo implantation of either the toric test lens (enVista one-piece hydrophobic acrylic toric IOL, Model MX60T) or the non-toric spherical control lens (enVista one-piece hydrophobic acrylic IOL, Model MX60) in one eye. Subjects in the higher astigmatic power cohorts (2.00 D, 2.75 D) were implanted with a test lens only in one eye. Postoperatively, subjects underwent complete ophthalmic evaluations at regularly scheduled intervals through Form 4 (Postoperative Days 120-180).

The test lens was the enVista toric IOL (Model MX60T). The effective corneal powers for each of the test lens plane cylindrical powers of the test IOLs are shown in **Table 1**.

Table 1:
enVista Toric IOL Cylinder Power

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¹Each Surgeon's individual surgically-induced astigmatism (SIA) was added to the recommended preoperative correction range to determine eligibility based on preoperative corneal cylinder. Once the SIA was estimated, this value stayed constant during the study for each investigator.

In order to facilitate toric IOL selection and axis placement, the Bausch + Lomb proprietary enVista toric Calculator was used to determine the appropriate enVista toric IOL model and axis of placement for each eye. The calculator accounted for surgically induced astigmatism

(SIA), incision location, and the subject's preoperative corneal astigmatism. In this trial all cataract incisions were to be placed on the preoperative keratometric steep axis.

Results

The results of the clinical study provide reasonable assurance that the Model MX60T IOL is safe and effective for the visual correction of aphakia and corneal astigmatism following cataract extraction.

The data support a significant dioptric reduction in cylinder and reduction in absolute cylinder, rotational stability of the lens, and improvement of both best corrected and uncorrected visual acuity at distance following implantation of the enVista toric IOL.

The primary effectiveness endpoints were mean toric IOL axial stability from Form 3 to Form 4, dioptric reduction in cylinder at Form 4, lens axis misalignment from surgical target markings at Form 4, and best corrected distance visual acuity at Form 4. All subjects in the toric IOL treatment groups demonstrated ≤5 degrees rotation from Form 3 (**Table 2**). Mean cylinder reduction from preoperative keratometric cylinder measurements in the randomized ITT population at Form 4 was 0.479 ± 0.665 D among those subjects with control IOLs and 0.865 ± 0.487 D among those subjects with 1.25 D toric IOLs (**Table 3**), showing a statistically significant improvement favoring the 1.25 D toric IOLs (P < 0.001). The mean percent reduction in absolute cylinder at Form 4 was 69.4 % for the all toric IOL Cohort and 36.8 % for the control IOL Cohort (**Table 4**). The percent of eyes within 0.50 D and 1.00 D of intended correction for All toric Cohort at Form 4 was 57.3 % and 90.9 %, respectively (**Table 5**). At Form 4, > 90 % of eyes in each toric IOL arm had misalignments of ≤ 10 degrees of intended markings, including 93.3 % of all toric IOL eyes (**Table 6**). Preservation of best-corrected distance visual acuity showed 99.1 % of eyes in the ITT population implanted with a toric IOL reported a VA of 20/40 or better at Form 4. Best-corrected distance visual acuity (BCDVA) results for the all toric IOL treatment group are presented in **Table 7** and **Table 8**. At Form 4, 109 subjects (99.1 %) in the All toric IOL Cohort achieved BCDVA of 20/40 or better. At Form 4, the mean ± SD UCDVA was 0.19 ± 0.16 logMAR in the control IOL treatment group and 0.11 ± 0.14 logMAR in the 1.25 D toric IOL treatment group (**Table 9**), which was a significant difference favoring the 1.25 D toric IOL arm (P < 0.001). At Form 4, 94.5 % of all toric IOL eyes and 83.3 % of control eyes had UCDVA of 20/40 or better.

The analysis of safety was based on the Safety cohort of 191 subject eyes for the implantation of a study lens (either test or control). The key safety outcomes are presented in **Table 10**. The rates of FDA defined potentially sight-threatening adverse events that occurred in the clinical trial at Form 4 were found to be less than the “FDA Grid” of historical controls. No serious adverse events occurred in the study eye. The incidence of adverse events experienced during the clinical trial was comparable to or lower than the incidence reported in the historic control (“FDA grid”) population (see **Table 10**). The enVista toric IOL demonstrated favorable safety compared with the Control IOL and the historical control (cf. ISO 11979-7 SPE) population, with no increase in incidence or severity of adverse events (AEs) compared with the Control IOL and no serious adverse events (SAEs) in the study eye. Overall, no safety signals were associated with the toric IOLs during this study.

Table 2:
Mean Toric IOL Axial Stability From Form 3 To Form 4 (ITT Population)

Absolute rotation (degrees)	Toric IOL			All Toric IOL (N=112)
	1.25 D (N=80)	2.00 D (N=20)	2.75 D (N=12)	
Absolute rotation from Form 3 to Form 4				
n	74	15	12	101
Mean ± SD	1.15 ± 1.08	0.92 ± 1.09	1.08 ± 0.73	1.11 ± 1.04
≤ 5 degrees rotation	74 (100.0 %)	15 (100.0 %)	12 (100.0 %)	101 (100.0 %)
95 % exact Confidence Interval	95.1 % to 100.0 %	78.2 % to 100.0 %	73.5 % to 100.0 %	96.4 % to 100.0 %

Table 3:
Mean Dioptric Cylinder Reduction From Preoperative Measurements (ITT Population)

Cylinder reduction (D)	Control IOL (N=79)	Toric IOL			All Toric IOL (N=112)
		1.25 D (N=80)	2.00 D (N=20)	2.75 D (N=12)	
Form 2					
n	77	80	20	12	112
Mean reduction ± SD	0.640 ± 0.591	0.966 ± 0.466	1.486 ± 0.498	2.115 ± 0.328	1.182 ± 0.594
Form 3					
n	79	79	18	12	109
Mean reduction ± SD	0.532 ± 0.627	0.864 ± 0.455	1.446 ± 0.519	1.926 ± 0.364	1.077 ± 0.584
Form 4					
n	78	80	18	12	110
Mean reduction ± SD	0.479 ± 0.665	0.865 ± 0.487	1.413 ± 0.532	1.944 ± 0.327	1.072 ± 0.601
Treatment effect at Form 4					
95 % Confidence Interval of effect		0.228 to 0.545			
Multiple imputation p-value ²		< 0.001			
² P-values and treatment effects are from a linear model Type II analysis, which include an effect for investigator and compare the control and 1.25 D toric IOLs at Form 4.					
Note: Dioptric change in cylinder = preoperative keratometric cylinder - postoperative manifest cylinder					

Table 4:
Mean Percent Reduction In Absolute Cylinder At Form 4

	Control IOL Mean ± SD	Toric IOL 1.25 D	Toric IOL 2.00 D	Toric IOL 2.75 D	All Toric IOL Mean ± SD
	N=79 n=78	N=80 n=80	N=20 n=18	N=12 n=12	N=112 n=110
Mean ³ % Reduction in Absolute Cylinder (± SD)	36.8 % ±50.49 %	64.8 % ±36.8 %	81.0 % ±31.3 %	82.8 % ±13.0 %	69.4 % ±34.8 %

³Mean of all subject (N) results at Form 4 - n (%)

Table 5:
Mean Percent Of Eyes With Reduction In Cylinder Within 0.50 D And 1.00 D Of Intended At Form 4 (ITT Population)

	Control IOL (N=79)	Toric IOL 1.25 D (N=80)	Toric IOL 2.00 D (N=20)	Toric IOL 2.75 D (N=12)	All Toric IOL (N=112)
Total Non-Missing, n	78	80	18	12	110
Within 0.50 D of Intended, n (%)	27 (34.6 %)	43 (53.8 %)	12 (66.7 %)	8 (66.7 %)	63 (57.3 %)
Within 1.00 D of Intended, n (%)	45 (57.7 %)	71 (88.8 %)	17 (94.4 %)	12 (100.0 %)	100 (90.9 %)

Table 6:
Mean Toric Lens Axis Misalignment From Surgical Markings At Form 4 (ITT Population)

Axis Misalignment	Toric IOL			All Toric IOL (N=112)
	1.25 D (N=80)	2.00 D (N=20)	2.75 D (N=12)	
Form 4 signed axis misalignment, degrees				
n	77	16	11	104
Mean ± SD	1.11 ± 8.69	3.08 ± 10.51	2.52 ± 3.03	1.56 ± 8.57
95 % tolerance interval for 90 % of the population	-15.51 to 17.73	-22.54 to 28.70	-5.78 to 10.82	-14.45 to 17.57
Form 4 absolute axis misalignment, degrees				
n	77	16	11	104
Mean ± SD	4.77 ± 7.33	5.15 ± 9.61	3.32 ± 2.01	4.68 ± 7.33
95 % tolerance interval for 90 % of the population	-9.25 to 18.79	-18.28 to 28.58	-2.19 to 8.83	-9.02 to 18.38
Form 4 absolute axis misalignment category, n (%)				
≤ 5 degrees	56 (72.7 %)	13 (81.3 %)	9 (81.8 %)	78 (75.0 %)
≤ 10 degrees	71 (92.2 %)	15 (93.8 %)	11 (100.0 %)	97 (93.3 %)
Form 4 signed axis misalignment category, n (%)				
-10.00 to -5.01 degrees	4 (5.2 %)	2 (12.5 %)	0 (0.0 %)	6 (5.8 %)
-5.00 to -0.01 degrees	31 (40.3 %)	3 (18.8 %)	2 (18.2 %)	36 (34.6 %)
0.00 degrees	0 (0.0 %)	0 (0.0 %)	0 (0.0 %)	0 (0.0 %)
+0.01 to +5.00 degrees	25 (32.5 %)	10 (62.5 %)	7 (63.6 %)	42 (40.4 %)
+5.01 to +10.00 degrees	11 (14.3 %)	0 (0.0 %)	2 (18.2 %)	13 (12.5 %)

Table 7:
Preservation Of BCDVA At Each Examination (All Toric IOLS, ITT Population)

BCDVA:	Preoperative	Form 2	Form 3	Form 4
	N=112	N=112	N=109	N=110
20/40 or Better, n (%)	27 (24.1 %)	108 (96.4 %)	109 (100 %)	109 (99.1 %)
Worse than 20/40, n (%)	85 (75.9 %)	4 (3.6 %)	0 (0.0 %)	1 (0.9 %)

Table 8:
BCDVA Without Glare At Form 4 (ITT Population)

	Control IOL (N=79)	Toric IOL 1.25 D (N=80)	Toric IOL 2.00 D (N=20)	Toric IOL 2.75 D (N=12)	All Toric IOL (N=112)
BCDVA (logMAR)					
Total Non-Missing, n	78	80	18	12	110
Mean (±SD)	0.01 (±0.09)	0.00 (±0.09)	0.05 (±0.10)	-0.01 (±0.09)	0.01 (±0.09)
BCDVA (Snellen)					
20/40 or Better, n (%)	78 (100.0 %)	79 (98.8 %)	18 (100.0 %)	12 (100 %)	109 (99.1 %)
Worse than 20/40, n (%)	0	1 (1.3 %)	0	0	1 (0.9 %)

Table 9:
UCDVA At Form 4 (ITT Population)

	Control IOL (N=79)	Toric IOL 1.25 D (N=80)	Toric IOL 2.00 D (N=20)	Toric IOL 2.75 D (N=12)	All Toric IOL (N=112)
UCDVA (logMAR)					
Total Non-Missing, n	78	80	18	12	110
Mean (±SD)	0.19 (±0.16)	0.11 (±0.14)	0.12 (±0.11)	0.13 (±0.18)	0.11 (±0.14)
UCDVA (Snellen)					
20/40 or Better, n (%)	65 (83.3 %)	76 (95.0 %)	18 (100.0 %)	10 (83.3 %)	104 (94.5 %)
Worse than 20/40, n (%)	13 (16.7 %)	4 (5.0 %)	0	2 (16.7 %)	6 (5.5 %)

Table 10:
ISO 11979-7 Safety And Performance Endpoints (SPE) Adverse Events Reported At Each Postoperative Visit, Implanted Subjects (Safety, All Toric IOL)

	Form 1 n/N (%)	Form 2 n/N (%)	Form 3 n/N (%)	Form 4 n/N (%)	Cumulative n/N (%)	ISO 11979-7 SPE rate (%)	p-value
Endophthalmitis	0/112	0/112	0/112	0/112	0/112	0.1	>0.999
Hypopyon	0/112	0/112	0/112	0/112	0/112	0.3	>0.999
Lens Dislocated From Posterior Chamber	0/112	0/112	0/112	0/112	0/112	0.1	>0.999
Cystoid Macular Edema	0/112	1/112 (0.9)	1/112 (0.9)	0/112	2/112 (1.8)	3.0	0.853
Pupillary Block	0/112	0/112	0/112	0/112	0/112	0.1	>0.999
Retinal Detachment	0/112	0/112	0/112	0/112	0/112	0.3	>0.999
Secondary Surgical Intervention	0/112	0/112	0/112	0/112	0/112	0.8	>0.999
Persistent Adverse Events: Noted at Form 4 ⁴							
Corneal Stromal Edema				0/112		0.3	>0.999
Iritis				0/112		0.3	>0.999
Cystoid Macular Edema				0/112		0.5	>0.999
Raised IOP Requiring Treatment				0/112		0.4	>0.999

⁴Cumulative versus persistent AEs are defined by the FDA SPE Grid and ISO 11979-7 as those occurring at Form 4 in this clinical study.

Adverse Event Reporting

Adverse events and/or potentially sight-threatening complications that may be regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be reported within five (5) days to Bausch & Lomb Incorporated. This information is being requested from all surgeons in order to document potential long-term effects of IOL implantation. Should any of these incidents take place in the European Union (EU), the event should also be reported to the competent authority of the EU member state in which the user is established.

Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with the IOL. These problems may be related to a specific lot of lenses or may be indicative of long-term effects associated with these lenses or with IOLs in general. If you wish to report a problem, in the USA, please call Bausch + Lomb at 1-800-338-2020. Outside

the USA, contact information can be found on www.bausch.com/contactus. In case of device explant/extraction from a patient, keep the device and contact Customer Service for return instructions.

How Supplied

Non-preloaded IOL:

The enVista IOL (model: ETE) is individually packaged in a sterile vial (containing blood bank saline) within a peel pouch, and should only be opened under sterile conditions.

Preloaded IOL:

The enVista IOL (models: ETPE, ETC20E, ETC24E) is preloaded in the SnapSet IOL shuttle and individually packaged in a sterile vial (containing blood bank saline), within a peel pouch, and should only be opened under sterile conditions.

A patient card and self-adhesive labels are supplied to provide traceability of the lens. The package is sterilized by gamma irradiation.

Expiry Date

Sterility is guaranteed unless the pouch is damaged or opened. The expiry date on the lens package is the sterility expiry date. This lens should not be implanted after the indicated sterility expiry date.

Safe Disposal

Dispose/discard of the unused or contaminated device/equipment and/or packaging by following applicable safe disposal procedures, and in accordance with applicable laws and regulations regarding the disposal of biohazardous materials.

Patient Registration Instructions And Reporting Registration

Each patient who receives an enVista IOL must be registered with Bausch + Lomb at the time of lens implantation. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to Bausch + Lomb. Patient registration is essential and will assist Bausch + Lomb in responding to adverse reaction reports and/or potentially sight-threatening complications. An implant identification card is supplied in the lens package and must be given to the patient.

Symbols And Abbreviations Used On Labeling

Symbol or Abbreviation	Symbol or Abbreviation Title
SE	Spherical Equivalent
CYL	Cylinder
IOL	Intraocular Lens
PC	Posterior Chamber
UV	Ultraviolet
D	Diopter
Ø _B	Body Diameter (Optic Diameter)
Ø _T	Overall Diameter (Overall Length)
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Fee paid for waste management
	Toric Monofocal Intraocular Lens


 Bausch & Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA

Manufactured by:
Bausch & Lomb Incorporated
21 Park Place Blvd. North
Clearwater, FL 33759 USA

®/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.
© 2025 Bausch & Lomb Incorporated or its affiliates

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