BAUSCH+LOMB

enVista ASPIRE™

HYDROPHOBIC ACRYLIC IOL

en

Device Description

The enVista Aspire™ hydrophobic acrylic intraocular lens (IOL) (non-preloaded model: EA / preloaded into shuttle models: EPA, EC20A, EC24A) was developed to replace the natural crystalline lens in adult patients in whom the cataractous lens has been removed. The composition and characteristics of the IOL are specified in the table below.

All IOL optical designs are associated with a certain amount of depth of focus. The Aspire IOL is classified as a monofocal IOL. Monofocal IOLs provide a limited depth of focus. The Aspire IOL uses an optical modification of the posterior aspheric surface to create a small continuous increase in IOL power within the central 1.5 mm diameter to slightly extend the depth of focus. However, clinically meaningful extension of the depth of focus has not been demonstrated in clinical trials.


In general, extending the depth of focus negatively affects the quality of vision at far distances. For the Aspire IOL, the power profile decreases towards the periphery outside the central 1.5 mm optic diameter, resulting in comparable distance image quality to the enVista monofocal IOL MX60E for a 3 mm pupil. Image quality has been estimated using non-clinical bench testing in an ISO2:+0.28 µm SA cornea model eye that simulates the spherical aberration from the cornea (see Figures 2 and 3).

The optic is designed with the SureEdge™ posterior-squared step edge to provide a 360-degree PCO barrier. The IOL employs an Accuset™ haptic with a broad, modified C-loop design and optic-haptic offset to facilitate improved contact and stability within the capsular bag. The enVista material that makes up the TruSight™ optic has been assessed for glistening-free capacity and scratch resistance. Incorporated StableFlex™ technology allows for enhanced IOL compliance for ease of loading, premium control through IOL delivery, and efficient optical recovery following implantation.

Physical Characteristics

Lens / Haptic Material	Hydrophobic acrylic (hydroxyethyl methacrylate (HEMA)-polyethylene glycol phenyl ether acrylate (poly(EG)PEA)-styrene copolymer, crosslinked with ethylene glycol dimethacrylate)
Material Characteristics	Index Of Refraction: 1.53 @ 35 °C; Specific Gravity: 1.19 g/ml
Optic Type / Powers	Aspheric / +6.0 to +34.0 Diopters in 0.5 Diopter increments (SE – Spherical Equivalent)
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 12.5 mm; Haptic Angle: 0°
Spectral Transmittance	**Ultraviolet**: UV (389) 10 % transmittance for +20.0 Diopter IOL. See Figure **1** with chart's **X** value = Wavelength (nm) and **Y** value = % Transmittance; chart compares the transmittance curve of an enVista IOL to a 53-Year-old human lens. **A** = +20 Diopter enVista EE IOL and **B** = 53-Year-Old Human Lens. **NOTE**: Light transmittance values for an IOL material may vary slightly depending on the method of measurement. **Reference**: 53-Year-old human lens data from Boettner, E.A. and Wolter, J. R., "Transmission of the Ocular Media," Investigative Ophthalmology, 1:776-783, 1962.

Figure 1

Indications

The enVista Aspire™ hydrophobic acrylic IOL (non-preloaded model: EA / preloaded into shuttle models: EPA, EC20A, EC24A) is indicated for primary implantation in the capsular bag of the eye in adult patients for the visual correction of aphakia following removal of a cataractous lens.

Warnings

As with any surgical procedure, there is risk involved. Physicians considering IOL implantation under any of the following circumstances should weigh the potential risk/benefit ratio:

1. Recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
4. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage to the endothelium during implantation.
6. Suspected microbial infection.
7. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.

Precautions For Use And Storage

1. Do not attempt to resterilize the IOL as this can produce undesirable side effects.
2. Prior to opening, inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use. The IOL should be used immediately after opening.
3. Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage (such as the loss of saline storage solution, or the presence of salt crystallization).
4. Store at room temperature. Do not freeze. Avoid high temperatures (>43 °C / >109 °F). Keep dry. Keep away from sunlight.
5. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution or sterile normal saline.
6. Do not place the IOL in contact with surfaces where such contamination can occur.
7. Do not autoclave the IOL.
8. Do not re-use the IOL. It is intended for permanent implantation. If explanted, sterility and proper function cannot be ensured.
9. Neither the safety and effectiveness, nor the effects of the Aspire IOL optical design on depth of focus, contrast sensitivity, and subjective visual disturbances (glare, halo, etc.) have been evaluated clinically. MTF testing of the Aspire IOL optical design (used in model EA) (Figures 2 and 3) may aid the surgeon in understanding the theoretical image quality expected with the Aspire IOL compared to the enVista monofocal IOL MX60E. However, these do not fully assess all aspects of clinical difficulties under all conditions. Surgeons must weigh the potential benefits of the modified optical design of the Aspire IOL (model EA) against the potential for risks associated with a degradation in vision quality and the lack of clinical data to characterize the impact of the Aspire IOL optical design on contrast sensitivity and subjective visual disturbance. These considerations may be especially relevant to patients with certain pre-existing ocular conditions (prior corneal refractive surgery, irregular corneal astigmatism, severe corneal dystrophy, macular disease, optic nerve atrophy, etc.) or intraoperative conditions (posterior capsular rupture, complications in which the IOL stability could be compromised, inability to place IOL in capsular bag, etc).
10. The safety and effectiveness of the IOL have not been substantiated in patients with pre-existing ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting an IOL in a patient with one or more of these conditions. Physicians considering IOL implantation in such patients should explore the use of alternative methods of aphakic correction and consider IOL implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.

 Before Surgery

 - Retinal conditions or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this IOL
 - Amblyopia
 - Clinically severe corneal dystrophy (e.g., Fuchs')
 - Rubella, congenital, traumatic or complicated cataracts
 - Extremely shallow anterior chamber, not due to swollen cataract
 - Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye (e.g., iritis or uveitis)
 - Aniridia
 - Iris neovascularization
 - Glaucoma (uncontrolled or controlled with medication)
 - Microphthalmos or macrophthalmos
 - Optic nerve atrophy
 - Previous corneal transplant
 - Pre-existing ocular conditions which may negatively impact stability of the implant

 During Surgery

 - Mechanical or surgical manipulation required to enlarge the pupil
 - Vitreous loss (significant)
 - Anterior chamber bleeding (significant)
 - Uncontrollable positive intraocular pressure
 - Complications in which the IOL stability could be compromised
11. Patients with preoperative problems such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, and chronic drug miosis may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from IOL implantation when such conditions exist.
12. A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on IOL implantation before attempting to implant the IOL.
13. As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block,

cyclitic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, acute corneal decompensation, toxic anterior segment syndrome (TASS), and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: IOL repositioning, lens replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.

14. Care should be taken to remove viscoelastic from the eye at the close of surgery.

Medical Device Re-Use Statement

If this product is reprocessed and/or re-used, Bausch + Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as 'single-use' which is defined as a device intended to be used once only for a single patient.

Calculation Of IOL Power

Suggested A-Constant: 119.1 (OPTICAL BIOMETRY)

The recommended A-Constant is intended for use with axial length measurements obtained by optical biometry. Use of axial length measurements by other techniques (e.g., Applanation A-scan) will normally require a different IOL constant. This number is a guideline only and is based on an evaluation of clinical data obtained using the IOL Master. The physician should determine preoperatively the power of the IOL to be implanted.

Directions For Use

1. Inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use.
2. Prior to implanting, examine the IOL package for type, power, and proper configuration.
3. Open the peel pouch and remove the vial in a sterile environment.
4. Remove the lid from the vial.
5. Follow steps below.
 a) Non-preloaded IOL (model: EA)
 1. With a pair of smooth forceps, remove the IOL from the vial by gently grasping the IOL haptic.
 2. Rinse the entire IOL with sterile balanced salt solution or sterile normal saline.
 3. Examine the IOL thoroughly to ensure particles have not become attached to it, and examine the IOL optical surfaces for other defects.
 4. The IOL may be soaked in sterile balanced salt solution until ready for implantation.
 5. It is recommended to use an approved inserter per the Validated Inserters table below.
 b) Preloaded IOL (models: EPA, EC20A, EC24A)
 1. Remove the SnapSet™ IOL shuttle from the vial by grasping the tab on the top.
 2. The SnapSet IOL shuttle may be flushed with sterile balanced salt solution during preparation for delivery.
 3. The SnapSet IOL shuttle preloaded with the enVista IOL (models: EPA, EC20A, EC24A) is designed for use with the EyeGility™ inserter for enVista preloaded.
6. It is recommended to use an approved viscoelastic for lubrication of the IOL during implantation. See table below.
7. There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.

Validated Inserters

Model	Inserter	Viscoelastic
EA	BLIS (BLIS-X1 cartridge)	Amvisc™ Plus, Amvisc™, OcuCoat™
	INJ100	
EPA (IOL only)	EyeGility (2.0 mm or 2.4 mm)	
EC20A (IOL + inserter)	EyeGility 2.0 mm	
EC24A (IOL + inserter)	EyeGility 2.4 mm	

Overview Of Clinical Studies

The enVista Aspire IOL model EA is a result of minor modifications (that did not require clinical study data) from the parent enVista model MX60. Clinical studies have not been conducted with the enVista Aspire IOL to assess the effect of its posterior aspheric surface on visual acuities, contrast sensitivity, or visual symptoms. Clinical studies have been conducted on the enVista single-piece IOL (model MX60) and the parent xact X-60 three-piece IOL (model X-60)*. The results of these studies are described herein.

* The AVS xact X-60 Intraocular Lens is not licensed for sale in Canada.

1. Summary Of Clinical Study For enVista Model MX60

A clinical study of the enVista Hydrophobic Acrylic Intraocular Lens, Model MX60, began in the United States on October 19, 2010. This prospective, single arm, open label study included a total of 122 subjects (122 eyes) at 6 clinical sites. Postoperatively, subjects underwent complete ophthalmic evaluations at regularly scheduled intervals through Form 4 (Postoperative Days 120-180).

Table 1 displays demographic information of subjects enrolled in the clinical trial. Table 2 displays BCVA results for best case subjects (those without clinically significant preoperative pathologies or macular degeneration at any time during the study) for 3 visits. At the Form 4 visit, 118 subjects (100 %) achieved BCVA of 20/40 or better, which exceeds the FDA grid of 96.7 %.

The key safety outcomes for this study are presented in Table 3. The rates of FDA defined potentially sight-threatening adverse events that occurred in the clinical trial at Form 4 were found to be less than the "FDA Grid" of Historical Controls. Two cumulative adverse events (2/122; 1.6 %) of cystoid macular edema were reported through the Form 4 visit. One persistent adverse event (1/121; 0.8 %) of cystoid macular edema was reported at the Form 4 visit. No serious ocular adverse events occurred during this study. One serious non-ocular adverse event of advanced leukemia with an outcome of death was reported during this study. The adverse event was determined by the study investigator to be unrelated to the investigational device, Model MX60 IOL. The incidence of adverse events experienced during the clinical trial was comparable to or lower than the incidence reported in the historic control ("FDA grid") population (see **Table 6**). The enVista IOL demonstrated favorable safety compared with the Control IOL and the historical control (cf. ISO 11979-7 SPE) population, with no increase in incidence or severity of adverse events (AEs) compared with the Control IOL and no serious adverse events (SAEs) in the study eye. Overall, no safety signals were associated with the IOLs during this study.

The results of clinical investigations provide reasonable assurance that the Model MX60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

Clinical Tables

Table 1:
Subject Demographics

	n	%
Number of Subjects	122	100.0
Gender		
Male	53	43.44
Female	69	56.56
Race		
African-American	1	0.82
Caucasian	119	97.54
Hispanic	2	1.64
Age		
< 60	13	10.66
60 to < 70	50	40.98
70 to < 80	54	44.26
≥ 80	5	4.10
Mean ± SD	69.0 (8.0)	
Range (Min, Max)	46, 93	

Table 2:
Best Corrected Visual Acuity By Postoperative Visit (Best Case Analysis Set)

Visual Acuity	Form 2		Form 3		Form 4	
	n	%	n	%	n	%
20/20 or better	82	69.5	93	78.2	99	83.9
20/25 or better	108	91.5	113	95.0	110	93.2
20/32 or better	117	99.2	117	98.3	117	99.2
20/40 or better	118	100	118	99.2	118	100
Exact 95 % CI for 20/40 or better[1]	96.9 %, 100 %		95.4 %, 100 %		96.9 %, 100 %	
Binomial Test p-value for comparing % of 20/40 or better, using FDA recommended p=0.967	1.0000					
20/41 to 20/63[2]	0	0.0	1	0.8	0	0.0
20/64 to 20/100	0	0.0	0	0.0	0	0.0
20/101 to 20/200	0	0.0	0	0.0	0	0.0
Worse than 20/200	0	0.0	0	0.0	0	0.0
N[3]	118		119		118	
N missing[4]	1		0		1	
Total N[5]	119		119		119	

[1]Exact 95 % confidence interval for the true percent from the study.
[2]Subject No. 04-027 experienced an AE of cystoid macular edema at Form 3; however, this AE was not persistent, and the subject achieved a BCVA of 20/20 at the Form 4 visit.
[3]Number of available best case subjects, plus best case subjects who missed scheduled visit but were seen later, and with VA measurement for the corresponding group [% = (n ÷ N) x 100 %].
[4]Number of subjects with missing VA measurement for the corresponding group.
[5]Total N: N + N-missing subjects for the corresponding group.

Table 3:
ISO Defined Cumulative And Persistent Adverse Events Through Form 4 (Safety Analysis Set)

Adverse Event	n/N[1]	%	ISO Grid (%)	90 % CI[2]	p-value[3]
Cumulative Safety Events[4]					
Endophthalmitis	0/122	0.0	0.1	0.00, 2.43	1.0000
Hyphema	0/122	0.0	2.2	0.00, 2.43	1.0000
Hypopyon	0/122	0.0	0.3	0.00, 2.43	1.0000
IOL Dislocation	0/122	0.0	0.1	0.00, 2.43	1.0000
Cystoid Macular Edema	2/122	1.6	3.0	0.29, 5.07	0.8839
Pupillary Block	0/122	0.0	0.1	0.00, 2.43	1.0000
Retinal Detachment	0/122	0.0	0.3	0.00, 2.43	1.0000
Secondary Surgical Intervention	0/122	0.0	0.8	0.00, 2.43	1.0000
Persistent Safety Events[5]					
Corneal Edema	0/121	0.0	0.3	0.00, 2.45	1.0000
Iritis	0/121	0.0	0.3	0.00, 2.45	1.0000
Cystoid Macular Edema	1/121	0.8	0.5	0.04, 3.86	0.4548
Raised IOP Requiring Treatment	0/121	0.0	0.4	0.00, 2.45	1.0000

[1]n: number of eyes reported with corresponding event. For cumulative event, N: number of implanted eyes. For persistent event, N: number of eyes returned for the Form 4 examination with non-missing response for the corresponding adverse event. A subject could be reported with more than one AE.
[2]Based on binomial distribution.
[3]Binomial test for the null hypothesis H_0: Percent from study ≤Percent from ISO Grid [per ISO 11979-7:2006 (E)].
[4]Occurring at any time during the study.
[5]Present at Form 4.

Other Clinical Findings

All subjects in the safety analysis set were evaluated for IOL glistenings at Form 3 and Form 4 visits. IOL glistenings were evaluated via retroillumination slit lamp examination utilizing a photographic grading scale provided in the protocol. The grading scale consisted of (in order of severity), "none, grade 0 (trace), grade 1, 2, 3, or 4." No glistenings of any grade were reported for any subject at any visit in the clinical study.

2. Summary Of Clinical Study For xact Model X-60 (Three-Piece IOL)

A clinical study of the xact Model X-60 IOL began in the United States on May 8th, 2002 and was conducted by Advanced Vision Science. A total of 383 subjects were enrolled, and 367 subjects were available for examination at one year, 312 were available at two years, and 281 were available at three years.

Table 4 displays demographic information of subjects enrolled in the clinical trial. Table 5 summarizes the best-corrected distance visual acuity (BCVA) results for best case subjects (those without clinically significant preoperative pathologies or macular degeneration at any time during the clinical trial).

Potentially sight threatening adverse events are listed in Table 6, along with the rate of occurrence in the clinical trial of the X-60 IOL, and are compared to the FDA Grid of Historical Controls. The number of patients included in the analysis of both cumulative and persistent adverse events in some cases was less than the number of patients who returned for examination and were available for analysis as a result of missing information in certain fields on the case report forms.

The results of clinical investigations provide reasonable assurance that the X-60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

Clinical Tables

Table 4:
Subject Demographics

	n	%
Number of Subjects	383	100.0
Gender		
Male	152	39.7
Female	231	60.3
Race		
Black	8	2.1
Caucasian	373	97.4
Hispanic	2	0.5
Age		
< 60	43	11.2
60 to < 70	105	27.2
70 to < 80	177	46.2
≥ 80	58	15.4
Mean ± SD	71.0 (9.11)	
Range (Min, Max)	45, 93	

Table 5:
Visual Acuity In Best Case Population

Visual Acuity	1 Year		2 Years		3 Years	
	n	%	n	%	n	%
20/20 or better	209	65.3	163	60.8	167	72.2
20/25 or better	275	85.9	215	80.2	203	86.3
20/30 or better	307	95.9	239	89.2	221	92.7
20/40 or better	317	99.1	253	94.4	229	95.1
FDA Grid for % of 20/40 or better	96.7 %		N/A		N/A	
N	320		268		242	

Table 6:
Cumulative And Persistent Adverse Events

Adverse Events	1 Year		FDA Grid 1 Year	2 Years		3 Years	
Cumulative Safety Events **Number of Eyes with Postop Visits = 382**	**n/N**	**%**	**%**	**n/N**	**%**	**n/N**	**%**
Endophthalmitis	0/382	0.0	0.1	0/382	0.0	0/382	0.0
Hyphema	0/382	0.0	2.2	0/382	0.0	0/382	0.0
Hypopyon	0/382	0.0	0.3	0/382	0.0	0/382	0.0
IOL Dislocation	0/382	0.0	0.1	0/382	0.0	0/382	0.0
Cystoid Macular Edema	3/376	0.8	3.0	3/377	0.8	3/377	0.8
Pupillary Block	0/382	0.0	0.1	0/382	0.0	0/382	0.0
Retinal Detachment	3/376	0.8	0.3	4/377	1.1	4/377	1.1
Secondary Surgical Intervention	1/382[1]	0.3	0.8	1/382	0.3	3/382[2,3]	0.8
Persistent Safety Events	**n/N**	**%**	**%**	**n/N**	**%**	**n/N**	**%**
Number of Eyes Available at the Visit	**367**		**N/A**	**312**		**281**	
Corneal Edema	0/366	0.0	0.3	0/312	0.0	0/281	0.0
Iritis	1/366	0.3	0.3	0/312	0.0	0/281	0.0
Cystoid Macular Edema	0/364	0.0	0.5	0/309	0.0	0/280	0.0
Raised IOP Requiring Treatment	0/366	0.0	0.4	0/312	0.0	0/281	0.0

[1]IOL was exchanged due to patient complaint of blurred vision, despite good BCVA. Investigator suspected glistenings might be related, however only modest improvement of vision was achieved after IOL exchange.
[2]IOL with glistenings was exchanged during retinal surgery to improve fundus visualization by the surgeon. Loss of vision was the result of retinal pathology and was not associated with the IOL.
[3]IOL was exchanged due to patient complaint of blurred vision. Investigator suspected glistenings might be related, however vision did not improve after IOL exchange. Since vision did improve after subsequent Nd:Yag capsulotomy, the complaint of blurred vision was not associated with the IOL.

Other Clinical Findings

In the IDE clinical trial, "glistenings" were observed in some cases. Glistenings, known to sometimes occur in some other hydrophobic acrylic IOLs, are microscopic vacuoles within the optic of the IOL that are visible through the slit lamp as multiple small refractile specks. Analysis of the clinical data confirmed no effect of glistenings on visual outcomes.

Testing established that glistenings were eliminated by a change in the IOL hydration solution from 10.0 % saline to 0.9 % saline. This was confirmed in an additional clinical trial conducted outside of the United States. In this study, 172 eyes of 142 patients were examined at least once between 1 and 6 months, and 123 eyes of 101 patients were examined at least once between 6 months and 2 years. No glistenings were observed at any time.

Adverse Event Reporting

Adverse events and/or potentially sight-threatening complications that may be regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be reported within five (5) days to Bausch & Lomb Incorporated. This information is being requested from all surgeons in order to document potential long-term effects of IOL implantation. Should any of these incidents take place in the European Union (EU), the event should also be reported to the competent authority of the EU member state in which the user is established.

Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with IOL. These problems may be related to a specific lot of lenses or may be indicative of long-term effects associated with these lenses or with IOLs in general. If you wish to report a problem, in the USA, please call Bausch + Lomb at 1-800-338-2020. Outside the USA, contact information can be found on www.bausch.com/contactus. In case of device explant/extraction from a patient, keep the device and contact Customer Service for return instructions.

How Supplied

Non-preloaded IOL:

The enVista IOL (model: EA) is individually packaged in a sterile vial (containing blood bank saline) within a peel pouch, and should only be opened under sterile conditions.

Preloaded IOL:

The enVista IOL (models: EPA, EC20A, EC24A) is preloaded in the SnapSet IOL shuttle and individually packaged in a sterile vial (containing blood bank saline), within a peel pouch, and should only be opened under sterile conditions.

A patient card and self-adhesive labels are supplied to provide traceability of the lens. The package is sterilized by gamma irradiation.

Expiry Date

Sterility is guaranteed unless the pouch is damaged or opened. The expiry date on the lens package is the sterility expiry date. This lens should not be implanted after the indicated sterility expiry date.

Safe Disposal

Dispose/discard of the unused or contaminated device/equipment and/or packaging by following applicable safe disposal procedures, and in accordance with applicable laws and regulations regarding the disposal of biohazardous materials.

Patient Registration Instructions And Reporting Registration

Each patient who receives an enVista IOL must be registered with Bausch + Lomb at the time of lens implantation. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to Bausch + Lomb. Patient registration is essential and will assist Bausch + Lomb in responding to adverse reaction reports and/or potentially sight-threatening complications. An implant identification card is supplied in the lens package and must be given to the patient.

Symbols And Abbreviations Used On Labeling

Symbol or Abbreviation	Symbol or Abbreviation Title
SE	Spherical Equivalent
CYL	Cylinder
IOL	Intraocular Lens
PC	Posterior Chamber
UV	Ultraviolet
D	Diopter
$Ø_B$	Body Diameter (Optic Diameter)
$Ø_T$	Overall Diameter (Overall Length)
℞ ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
	Fee paid for waste management
	Monofocal Intraocular Lens

Bausch & Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA

Manufactured by:
Bausch & Lomb Incorporated
21 Park Place Blvd. North
Clearwater, FL 33759 USA

®/TM are trademarks of Bausch & Lomb Incorporated or its affiliates.
© 2025 Bausch & Lomb Incorporated or its affiliates

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